GREINER
GLASINSTRUMENTE
GMBH

Refraktometer GREINORM

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
GREINORM Brix Refractometer

Operation Instruction

1. Prism 2. Cover plate 3. Calibration 4. Mirror tube 5. Eyepiece (Adjusting ring of diopter)

Method of operation

Point the front of the refractometer in the direction of a bright source of light and adjust it with help of the adjusting ring of diopter (5) until the reticle can be seen clearly.

Calibration

The refractometers are delivered calibrated. New calibration is necessary if the instrument got a shock (by transport or similar) or if it is not used a longer time. In this case please open the cover plate of the instrument and drop one or two drops of distilled water on the prism. Then close the cover plate and press it gentle. Afterwards adjust the correct screw (3) to make the light/blue boundary coincide with zero line.

Attention

- ☞ Adjustment for zero liquid and specimen should be in the same temperature
- ☞ Do not use water to wash the instrument
- ☞ Do not touch or scratch the optical surface and avoid strong shocks

Temperature Compensation

The reference temperature is 20°C. ATC is mounted with Automatic Temperature Compensation, it is the extra function which enables user concentrate on measurements without worrying the temperature. 10 till 30°C is the possible range of compensation.

Care

You should clean the measuring prism before and after using carefully with a clammy and soft cloth and after that it should be rubbed dry. The prism consists of a high-class glass and therefore it is very sensible. If the testing liquid will be left on the glass without cleaning it can attack the glass. The same valids also for water. The prism and the coverplate should be cleaned with a soft cloth accurate and the prism should be dry after this else it will also attack the prism and coverplate. **The prism can be cleaned with pure alcohol.** The measuring prism may not be scratched. The instrument must be stored dry and dust-free. As expensive cameras the refractometers do not suffer rude handling. You should save it for falling down and kicks.

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Greiner Glasinstrumente GmbH